



Pharmaceutical Test List

10329 Stony Run Lane - Ashland, VA 23005
 800-888-8061 - 804-365-3002 (Fax)
www.analyticscorp.com

Compound	Media	ng/ Sample	Method
5-Fluorouracil	J	10	HPLC-UV
Acetaminophen	E, J	1000	HPLC-UV
Acetylsalicylic acid	E, J	2000	HPLC-UV
Aciclovir	E, J	9000	HPLC-UV
Acitretin	J	50	HPLC-UV
Actinomycin D	Wipe	100	HPLC-UV
Amitriptylene	E, J, A	100	HPLC-UV
Amoxicillin	E, J	100	HPLC-UV
Amphotericin B	J	400	HPLC-UV
Anagrelide Hydrochloride	N	70	HPLC-UV
Atracurium besylate	E, J	600	HPLC-UV
Atropine Sulfate	H	125	HPLC-UV
Azothioprine	I/Wipe	100/1000	HPLC-UV
Bendamustine	H/Wipe	100/1000	HPLC-UV
Benzoate	A	10	HPLC-UV
Bromocriptine mesylate	E	500	HPLC-UV
Bupropion HCl	E, J	400	HPLC-UV
Butalbital	E	2000	HPLC-UV
Butoconazole Nitrate	A	500	HPLC-UV
Capsaicin	J	20	HPLC-UV
Capecitabine	E, J	100	HPLC-UV
Cefaclor	Wipe	200	HPLC-UV
Cefamandole Nafate	Wipe	100	HPLC-UV
Cefazolin Sodium	Wipe	100	HPLC-UV
Ceftazidime	Wipe	100	HPLC-UV
Ceftriaxone sodium	N/A	N/A	HPLC-UV
Cefuroxime Axetil	Wipe	100	HPLC-UV
Cefuroxime Sodium	Wipe	200	HPLC-UV
Cephalexin	Wipe	100	HPLC-UV
Cephalothin	Wipe	300	HPLC-UV
Chlordiazepoxide	E, J	40	HPLC-UV
Chlormadinone Acetate	A	50	HPLC-UV
Chlorohydrine	A	1000	HPLC-UV
Clenbuterol	E, J	20	HPLC-UV
Clidinium bromide	E, J	2000	HPLC-UV
Clobetasol propionate	E, J	100	HPLC-UV
Clonazepam	E, J	200	HPLC-UV
Clozapine base	E	780	HPLC-UV
Cyclophosphamide	E, J, A	8000	HPLC-UV

Compound	Media	ng/ Sample	Method
Cyclosporin A	E	500	HPLC-UV
Darifenacin	F	125	HPLC-UV
Desogestrel	E	20	HPLC-UV
Diatrizoic Acid	E	156	HPLC-UV
Digoxin	N/Wipe	150	HPLC-UV
Dihydrocapsaicin	J	20	HPLC-UV
Doxifl Uridine	J	25	HPLC-UV
Doxyclyne Hyclate	F	500	HPLC-UV
Ethinyl Estradiol	E/Wipe	10/60	HPLC-FL
Etretinate	K	50	HPLC-UV
Floxuridine	J	50	HPLC-UV
Flucytosine	J	40	HPLC-UV
Fludrocortisone Acetate	H	10	HPLC-UV
Flumethasone	A	10	HPLC-UV
Fluocinolone Acetonide	A	20	HPLC-UV
Fluocinonide	A	20	HPLC-UV
Fluorescein	L	625	OD
Fluorouracil (5-)	E, J	100	HPLC-UV
Flurazepam	E, J	200	HPLC-UV
Fluticasone Propionate	E, J	2000	HPLC-UV
Ganciclovir Sodium	A	100	HPLC-UV
Hydrazoic Acid	M	50 (300)	HPLC-UV
Hydrocodone	J	200	HPLC-UV
Hydrocodone Bitartrate	J	200	HPLC-UV
Insulin	F	2000	HPLC-UV
Isotretinoin	J	100	HPLC-UV
Ketorolac Tromethamine	A	100	HPLC-UV
Labetalol	E, J	1000	HPLC-UV
Lactose	A/ Wipe	2	IC
Laidlomycin	A	400	HPLC-UV
Laidlomycin Propionate	A	400	HPLC-UV
Lamivudine	E, J	1000	HPLC-UV
Lansoprazole	E	30	HPLC-UV
Lasalocid	E, J	800	HPLC-UV
Levonorgestrel	E	10	HPLC-UV
Mannitol	A/Wipe	1	IC
Mestranol	E, J	80	HPLC-UV
Methamphetamine	E, J	200	HPLC-UV
Methylphenidate	E, J	200	HPLC-UV

Compound	Media	ng/ Sample	Method
Methylphenidate (dex-)	E, J	200	HPLC-UV
Mivacurium HCl	E, J	800	HPLC-UV
Morphine Sulfate	E, J	100	HPLC-UV
Mycophenolate Mofetil	A	200	HPLC-UV
Nadolol	E, J	10	HPLC-UV
Naproxen	A	2	HPLC-FL
Naproxen Sodium	A/Wipe	2	HPLC-FL
Nicardipine Hydrochloride	A	1000	HPLC-UV
Norethindrone	A	10	HPLC-UV
Norethindrone Acetate	E/Wipe	30	HPLC-UV
Norgestimate	E	10	HPLC-UV
Orlistat	B	1000	HPLC-UV
Oseltamivir	E, J	100	HPLC-UV
Oxfendazole	A	1000	HPLC-UV
Oxytetracycline	F	2500	HPLC-UV
Pergolide Mesylate	E	8.5	HPLC-FL
Pilocarpine HCl	E	600	HPLC-UV
Prasterone (DHEA, dehydroepiandrosterone)	E	100	HPLC-UV
Procainamide HCl	E, J	100	HPLC-UV
Procarbazine Hydrochloride	J	100	HPLC-UV
Progesterone	A	10	HPLC-UV
Pyridostigmine Bromide	E, J	200	HPLC-UV
Ramoxifen Citrate	G	200	HPLC-UV
Ranitidine HCl	J	400	HPLC-UV
Ranolazine Free Base	A	500	HPLC-UV
Ranolazine Hydrochloride	A	500	HPLC-UV
Retinol	O	100	HPLC-UV
Ribavirin	H	30	HPLC-UV
Rifaximin	E	17,000	HPLC-UV
Salinomycin Sodium	J	500	HPLC-UV
Sertindole	E	200	HPLC-UV
Sumatriptan Succinate	E, J	1000	HPLC-UV
Temazepam	E	500	HPLC-UV
Testosterone Propionate	A	10	HPLC-UV
Tezacytabine (FMdC)	J	20	HPLC-UV
Thalidomide	C	30	HPLC-UV
Ticlopidine Hydrochloride	A	200	HPLC-UV
Tobramycin	E	2000	HPLC-UV



Pharmaceutical Test List

10329 Stony Run Lane - Ashland, VA 23005
800-888-8061 - 804-365-3002 (Fax)
www.analyticcorp.com

Compound	Media	ng/ Sample	Method
Tolcapone	B	5000	HPLC-UV
Trenbolone Acetate	A	5	HPLC-UV
Tretinoin, Vesanoid	J	50	HPLC-UV
Triamcinolone	K	50	HPLC-UV
Valaciclovir	J	3000	HPLC-UV

Compound	Media	ng/ Sample	Method
Verapamil	J	200	HPLC-UV
Vitamin D3	D	120	HPLC-UV
Warfarin	E	50	HPLC-UV
Zalcitabine	A	100	HPLC-UV
Zidovudine	E, J	3000	HPLC-UV

Compound	Media	ng/ Sample	Method
Zopiclone	E	38	HPLC-UV

Media Codes:

- A** - 37 mm glass filter type A/E (1µm) and pad in polystyrene cassette
- B** - 37 mm quartz filter (1µm) and pad in polystyrene cassette
- C** - Two 37 mm quartz filters (1µm) in polystyrene cassette
- D** - 37 mm Teflon filter (1µm) and pad in polycarbonate cassette
- E** - 37 mm Teflon filter (1µm) and pad in polypropylene cassette
- F** - 37 mm Teflon filter (1µm) and pad in 3-piece polypropylene cassette
- G** - 37 mm Teflon filter (1µm) and glass filter type A/E as support pad in polypropylene cassette
- H** - 25 mm Teflon filter (1µm) and pad in polypropylene cassette
- I** - 37 mm Teflon filter (0.45µm) and pad in polypropylene cassette
- J** - 37 mm Teflon filter (1µm) and pad in polystyrene cassette
- K** - 37 mm Teflon filter (0.5µm) and pad in polystyrene cassette
- L** - Two 37 mm quartz filters (1µm) in polypropylene cassette
- N** - Two 37 mm glass filter type A/E (1µm) in polystyrene cassette
- O** - 25 mm impregnated Teflon filter (1µm) in IOM cassette