

HEMODIALYSIS LABORATORY TEST REQUEST

ACCOUNT NAME AND ADDRESS



10329 Stony Run Lane Ashland, VA 23005
 Phone (804) 365-3000
 Toll Free (800) 888-8061
 Fax (804) 365-3002
 www.analyticscorp.com

DATE SHIPPED
MACHINE NUMBER (SAMPLE ID)
FACILITY NAME
SAMPLE DATE
TYPE OF WATER (CLASSIFICATION) DO NOT SEND DIALYSATE

CHECK ONE:

AAMI ANALYSIS
 ANSI/AAMI RD62-2006

Aluminum	Chromium	Potassium
Antimony	Copper	Selenium
Arsenic	Fluoride	Silver
Barium	Lead	Sodium
Beryllium	Magnesium	Sulfate
Cadmium	Mercury	Thallium
Calcium	Nitrate	Zinc

SPECIAL INSTRUCTIONS AND/OR UNUSUAL CONDITIONS:

LAL/Endotoxin (Special Kit)

HPC

TOC

Other

CONTACT NAME	TELEPHONE NUMBER	EMAIL/ FAX NUMBER
--------------	------------------	-------------------

HEMODIALYSIS WATER SAMPLE KIT INSTRUCTIONS

Each kit contains one bottle one zip lock plastic bag, and one pre-labeled mailer.

1. Run water tap or access line for two (2) minutes.
2. Rinse bottle with a small amount of water to be sampled.
3. Fill bottle completely with water sample, cap tightly, and place filled bottle in the plastic bag.
4. Fill out the Hemodialysis Laboratory Test Request form.
5. Place bottle enclosed in plastic bag and completed test request inside mailing container and mail.

*Ship LAL/HPC analysis cold but not frozen overnight, first morning delivery to preserve the integrity of the sample.
 If you have any questions, please call Customer Services at 800-888-8061.

PLEASE COMPLETE A SEPARATE FORM FOR EACH WATER SAMPLE